BIODAM should not be used in patients is known to have sensibility to any of the ingredients listed. should not be used in the following cases:

- Patients who already have gingivitis or periodontal illness and gingival injuries;
- Patients who had a periodontal surgical intervention recently;


WARNINGS:

Do not use material after expiry date.

The manufacture is not responsible for any damage resulting from a misuse or not provided for the instructions of use.

“KEEP MATERIAL OUT OF THE REACH OF CHILDREN”.

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Reg. ANVISA: 10298550101.

COMPOSITION:

Dimethacrylate groups; silicon dioxide; BHT; catalysts; inorganic filler and pigments.

INDICATION:

BIODAM is indicated as a gingival protector in cases of use products that offer risks of irritation (For ex: bleaching agents based on hydrogen peroxide for clinical use), preventing the contact of the irritating agent with the gingival tissue.

TECHNICAL INFORMATION:

For the professional bleaching procedure with hydrogen peroxide or for the gingival isolation from any other product, is necessary to use a material that offers protection. **BIODAM** is an insulator and protector material. After lighting cured creates a rigid mechanical barrier to isolate the teeth and protect the gingival tissues. **BIODAM** can substitute or complement the rubber dam for the gingival protection.

BIODAM promotes a relative isolation, adhering to the gingiva and mechanically blocking any product crossing. **BIODAM** is ease to apply and its green color facilitates the contrast with the teeth. Adheres well to the gingiva and its removal is easily done with the help of an instrument (For ex: explorers), detaching the product.

USE INSTRUCTIONS:

1. Dry well the: Teeth and the gingiva, especially the gingival portion;
2. Apply a thin layer of **BIODAM**, approximately from 4 to 6 mm of height and from 1 to 1,5 mm of thickness on the gingiva. Cover about 0,5 mm of the enamel, in the cervical region, to prevent the contact with the gingival border and repeat the process in all teeth that will be isolated.
Extensions of **BIODAM** between the teeth gives greater retention of protection.
3. Light cure for 20 to 30 seconds and continue this procedure until all length of **BIODAM** has been light cured.
4. Verify if the product is adequately light cured and if it offers protection on the all gingivas surface needed.
5. For the removal, detach the product onward with the help of an instrument.

SPECIAL CARES:

Light cure the product in accordance with the use instructions, respecting the time and the irradiation place. The incomplete or an extreme light cure may cause a light gingival irritation.

Dispose the product in accordance with local regulations, ensuring its complete distortion, preventing its reuse and environment damages

CONTRAINDICATION:

biodinâmica