
# VARNAL

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biodinâmica®

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Reg. ANVISA: 10298550045

### COMPOSITION:

Resin Staybilite (Ester 10) and Dimethylcetone.

### INDICATION:

Performs as a barrier against the chemical elements that compose the restoring material; Reduce marginal leakage around most of restoring material, specially in amalgam restorations. With such reduction in the marginal leakage, it acts in the control of inflammatory reaction and in the reduction of post-surgical sensitivity;

Used in the galvanic shock treatment.

### TECHNICAL INFORMATION:

VARNAL is a cavity varnish used as auxiliar in covering and protection for the pulp-dentin complex.

When VARNAL is applied there is a solvent evaporation an a thin resin film remains on the surface of cavity prepare. Such thin film acts as a semi permeable membrane, obstructing the passage of some ions. It is most important to obtain an uniform and continuous covering on the whole surface of cavity prepare. If the layer is uneven or has pores, the result will not be satisfactory.

### USE INSTRUCTIONS:

1. VARNAL must be applied aiming to obtain an uniform and continuous layer on all over surface of the prepared cavity. Such application should be done with a disposable brush or a cotton ball; NOTE: The disposable brush or cotton ball must not stay dribbling with varnish excess.

2. We recommend two successive applications, to reduce the possibility of porus formation and to provide a continuous layer. Notice that a 15 to 20 seconds interval is required between the applications, to let the first layer of varnish dry.

Whenever necessary, the product can be diluted with an appropriate solvent (dimethylcetone) in order to provide the desired film thickness.

### SPECIAL CARES:

Product must be kept inside hermetically sealed flasks, in fresh places sheltered from sunlight and away from fire. When in use keep it opened just the strictly necessary time until obtain the portion needed in order to avoid evaporation of its solvent.

Dispose the product in accordance with local regulations, ensuring its complete distortion, preventing its reuse and environment damages.

### CONTRAINDICATION:

VARNAL should not be used in patients is known to have sensibility to any of the ingredients listed.

VARNAL should not be used on acrylic restorations or composite restorations, because the conventional cavity varnishes interferes with polymerization reactions.

### WARNINGS:

Do not use material after expiry date.

The manufacturer is not liable for damages caused by improper use or use not provided at these instructions.  
"KEEP MATERIAL OUT OF THE REACH OF CHILDREN".

FOR USE IN DENTISTRY ONLY