ES

IT

DE

EN

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General Precautions:

• U.S. Federal law restricts this device to sale by or on the order of a dentist.
 • Review instructions and SDS (www.dentaltech.com/resource) for all provided materials prior to beginning treatment. Use only as directed.
 • Do not use in the presence of eugenol based materials as complete setting may not occur.
 • Lower temperatures will increase the setting time and viscosity.
 • Mixing tips and intra oral tips (if supplied) are disposable and intended for single patient use.

Adverse Reactions:

Allergic or pre-existing skin, eye and respiratory conditions may be aggravated by exposure to this product

Instructions for Use:

1. Isolate the tooth to be prepared by using either a rubber dam or cotton rolls.
2. Remove all existing restorations and decay from the tooth.
3. Apply **Alpha-dent®** etchant or **Alpha-Etch37®** to the occlusal surface following the etchant instructions. Wait for 15 seconds; never leave etchant on tooth for longer than 60 seconds. Lightly agitating the etchant during its application can enhance the action.

4. Rinse thoroughly with water. Do not dry the surface. Remove excess water with a brief burst of air, leaving the dentin and enamel visibly moist.
 5. Apply bonding agent according to manufacturer's instructions. We recommend the use of **Alpha-Bond®** Adhesive with Activator.
 6. Select the desired shade and place the cartridge into the Type 25ml Applicator gun. Review cartridge assembly instructions below.
 7. Dispense with slow and steady pressure, keeping the tip at the base of the preparation allowing the material to rise around the tip to eliminate voids. If desired for more accurate placement, an inter-oral tip can be inserted onto the open end of the mixing tip.
 8. **Alpha-Core® DC** has a 1 minute 30 second working time and 5 minute total setting time if used as a self cure material. If used as a light cure material: using a visible light-curing unit (minimum output of 600 mW/cm²) cure the material for 40 seconds. You can expect a maximum depth of cure to be 4mm. At this layer (<2mm) needs to be cured for only 20 seconds when using a visible light-curing unit (minimum output of 600 mW/cm²).
 9. Replace cap immediately after use by aligning the notch on the cap with the corresponding notch on the cartridge. Once aligned, push cap onto the cartridge and rotate the cap until it locks in position (90°).

10. Finish **Alpha-Core® DC** when fully set.

Cartridge Assembling

1. Make sure that the plunger slide is pulled back completely. This can be accomplished by pushing up the small lever below the plunger slide while pulling back the plunger slide.
2. Lift the plastic lever on the top of the gun and insert the cartridge, making sure that the notch on the bottom of the cartridge aligns with the notch on the applicator gun. Push down firmly on the lever to lock cartridge into place.
3. Remove the cartridge cap by rotating the cap 90° and then pulling.
4. To ensure an even flow, when using a new cartridge for the first time, extrude a "pea" size amount of material and discard. Be sure to avoid cross contamination of the pastes as setting within the cartridge may occur.

5. Place the mixing tip onto the cartridge by aligning the notch on the mixing tip with the corresponding notch on the cartridge. Once aligned, push the mixing tip onto the cartridge and rotate the mixing tip until it locks into position (90°).

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