TECHNICAL INFORMATION:

The bonding agent SILANO acts as a connector between the indirect restoration and the composite by means of a stable chemical bond, providing increased bond strength in the indirect/composite restoration interface.

USE INSTRUCTIONS:

FOR CEMENTATION OF INDIRECT RESTORATIONS WITH RESINOUS CEMENTS:

1. Prepare the tooth surface according to the technique recommended by the agent cement;
2. Apply porcelain etchant (hydrofluoric acid, e.g. Ataque F) to the inner surface of the restoration for the time recommended by the manufacturer. Wash and dry carefully;
3. Apply SILANO BONDING AGENT using a brush on the conditioned surface, wait 1 minute;
4. Cement the restoration according to the recommendations for the resin cement to be used (ex Master Cement Dual).

FOR COMPOSITE RESIN REPAIRS:

1. Prophylaxis of the surface to be repaired with pumice stone and water, dry and isolate;
2. Remove friable edges and roughen the surface with diamond burs;
3. Etch the surface with porcelain conditioner (hydrofluoric acid, e.x. Ataque F) at least time recommended by the manufacturer;
4. Wash thoroughly and dry completely;
5. Apply the SILANO JOINING AGENT using a brush over the entire surface conditioned, wait 1 minute;
6. Apply the selected conventional adhesive system (ex Master Bond) and light cure as required. manufacturer's instructions;
7. Repair the fractured area with a composite (eg Master Fill).

COLLAGE OF FRAGMENTS:

1. Perform prophylaxis with pumice stone and water on the surface to be repaired, dry and isolate;
2. Roughen the surface to be repaired and the fragment with diamond burs;
3. Condition the surface of the fragment and the flaw with porcelain conditioner (acid hydrofluoric, e.g. Attack F) for the time recommended by the manufacturer;
4. Wash thoroughly and dry thoroughly;
5. Apply the SILANO JOINING AGENT using a brush over the entire surface conditioned, wait 1 minute;
6. Cement the fragment according to the cement manufacturer's instructions;
7. If there is a lack of material in the fragment (ceramic, porcelain, resin), rebuild the

fractured area with composite (ex Master Fill).

SPECIAL CARES:


- Every procedure must be done preferably with absolute isolation;
- Extensive repairs with composites should be performed in layers to reduce the effect of polymerization contraction.
- Dispose of the product in accordance with local legislation, ensuring its complete mischaracterization in order to prevent its reuse and environmental damage.

CONTRAINDICATION:

Contraindicated for people with reported sensitivity to any of the product's components. Do not use material after expiry date. The manufacturer is not liable for damages caused by improper use or use not provided at these instructions.

"KEEP MATERIAL OUT OF THE REACH OF CHILDREN".
FOR USE IN DENTISTRY ONLY

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