

TECHNICAL INFORMATION

INTERIM Temporary Restorative Material

CONTENTS:

The product is presented on powder and liquid form; 1x38g INTERIM Powder; 1x15mL INTERIM Liquid.

COMPOSITION:

INTERIM POWDER: Zinc Oxide and Acrylic Resin; INTERIM LIQUID: Eugenol and Acetic Acid Glacial.

INDICATION:

INTERIM is a temporary intermediate restorative material; It is used as therapeutic basis for temporary restoration and cavities filling.

TECHNICAL INFORMATION:

Setting time is approximately 5 to 8 minutes (for powder / liquid 6/1 in weight). It's hydrogenic concentration has pH = 7 - 8, even at the moment that the product is placed in tooth. By this reason INTERIM is one of the less irritant cement used in Dentistry.

USE INSTRUCTIONS:

Powder/liquid ratio (g/g): 6/1 in weight.

1. On a flat glass plate dispense one measure of powder and one drop of liquid, divide powder into 3 equal parts:

2. Mix 1/3 of the powder with the liquid and espatule during 20 seconds;

3. Add more 1/3 of the powder and espatule during 20 seconds;

4. Finally add remaining 1/3 of the powder and espatule during 20 seconds until obtain a low adherent mass to the glass plate:

5. Place the material using Hollemback plugger and finish placement with a plugger or a burnisher while the cement is setting. Any excess may be taken off using a Hollemback carver or low speed burs.

CONTRAINDICATION:

INTERIM should not be used in patients is known to have sensibility to any of the ingredients listed.

CHARACTERISTICS/ PROPERTIES:

Zinc Oxide and Eugenol (modified) based material: Reduced setting time - 5 to 8 minutes: The color of the powder is ivory; Easy to manipulate and insertion inside dental cavity Presents strength:

Presents resistance to abrasion;

biodinâmica respect for you

TECHNICAL INFORMATION

Antiseptic and sedative properties due to Eugenol; Insoluble on oral environment; Biocompatibility as a pulp protective Promotes marginal sealing; Easy to remove.

BIBLIOGRAFHY:

- ✓ SKINNER, E.; MATERIAIS DENTÁRIOS de SKINNER: Tradução Dioracy Fonterrada Vieira;
- ✓ 8ªEd. Rio de Janeiro: Editora Interamericana, 1984;
- ✓ MONDELLI, J. e outros; DENTÍSTICA OPERATÓRIA 3ª Edição, São Paulo, Sarvier, 1977;
- ✓ BERBERT, A., 1934 Endodontia Prática São Paulo: Sarvier, 1980;
- ✓ VIEIRA, D.F.- BASES PARA A APLICAÇÃO RACIONAL DOS MATERIAIS ODONTOLÓGICOS. 2ª Ed. São Paulo, Atheneu, USP, 1976;
- SSAO, M.; & PINTO, A. C. G.; MANUAL DE ODONTOPEDIATRIA Editora Dent'Art S.A. São Paulo 1973.